

**ROBIC SC-500E
Operating Instructions**



To Time Single Events

Press “A” to **START** timing.

Press “A” to **STOP** timing.

After reviewing the time reading, press “B” to **RESET** timer to all zeros.

Note: You may take as many consecutive single event times as necessary by repeating this procedure.

Silent or Audible Operation

Pressing “C” will change from silent to audible operation. You may change this function at any time.

The SC-500E will time up to 10 hours with 1/100 second resolution.

PLEASE NOTE: The “D” (right front) pushbutton DOES NOT provide any function. Please avoid pressing this button.